COATEST® SP FVIII - 82 4086 63

Intended use

For the photometric determination of factor VIII activity in citrated plasma, such as when identifying factor VIII deficiency or monitoring patients on replacement therapy, as well as for potency estimation of FVIII concentrates.

Measurement principle

Composition
 1. **S-2765 15.4 mg + I-2581 1 vial**
 Chromogenic substrate (N-a-Z-D-Arg-Gly-Arg-pNA), 15.4 mg, synthetic thrombin inhibitor, 0.4 mg, and manniitol added as a bulking agent.
 2. **Factor IXa + factor X 9.2 IU 1 vial**
 Lyophilized bovine factors IXa and X with bovine albumin added as a stabilizing agent.
 3. **CaCl₂ 6 ml 1 vial**
 Calcium chloride solution, 0.025 mol/L
 4. **Buffer, stock solution 20 mL 1 vial**
 20 mL concentrated Tris buffer containing NaCl and BSA. Characteristics of tenfold diluted buffer: Tris 0.05 mol/L, pH 7.3, 10 mg/L Ciprofloxacin and 1.0% BSA.
 5. **Phospholipid 2 mL 1 vial**
 Mixture of highly purified phospholipids and 10 mg/L Ciprofloxacin.

PRECAUTION AND WARNINGS
Danger
FIXa + Factor X 9.2 IU Hazard class: Resp. sens. 1, H334
Hazard statements: H334: May cause allergy or asthma symptoms or breathing difficulties if inhaled.
Precautionary statements: P261: Avoid breathing vapors/ spray; P284: [In case of inadequate ventilation] wear respiratory protection; P304 + P340: IF INHALED: Remove person to fresh air and keep comfortable for breathing; P305 + P311: If experiencing respiratory symptoms: Call a POISON CENTER/doctor. P501: Dispose of contents/container in accordance with local/regional/national/international regulation.
Supplemental hazard information: Contains Factor IXa. Up to 2.95 % of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and for the aquatic environment.
CaCl₂ S-2765 15.4 mg + I-2581, Buffer , Phospholipid
Hazard class: none
Hazard statements: none
Precautionary statements: none
Supplemental hazard information:
CaCl₂ EUH210: Safety data sheet available on request.
S-2765 15.4 mg + I-2581: EUH210: Safety data sheet available on request. Up to 6.14% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and for the aquatic environment.
Buffer: EUH210: Safety data sheet available on request. Up to 7 % of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health.
Phospholipid: Up to 40.4 % of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health.
 This product is for *in vitro* diagnostic use.

Preparation
 The reagents are reconstituted according to the specific instrument application. For microplate and test tube techniques:
 1. **S-2765 + I-2581:** Reconstitute with 12.0 mL of sterile water or NCCLS type II water¹¹, to obtain a concentration of 1.8 mmol/L.
 2. **Factor IXa + Factor X:** Reconstitute with 10.0 mL of sterile water or NCCLS type II water¹¹
 3. **CaCl₂:** Ready to use.
 4. **Buffer, stock solution:** Dilute 1:10 (1+9) with sterile water or NCCLS type II water¹¹. Prepare a new buffer working solution each day.
 5. **Phospholipid:** Ready to use.

Reagent storage and stability
 When kept at 2-8°C the sealed reagents are stable until the expiry date printed on the label. Contamination by microorganisms should be avoided once the vials are opened.
 1. **S-2765 + I-2581:** Stability after reconstitution: 3 months at 2-8°C.
 2. **Factor IXa + Factor X:** Stability after reconstitution: 12 hours at 2-8°C. The solution can be stored frozen in aliquots at -20°C (or at lower temperature) for 3 months. Do not refreeze.
 3. **CaCl₂ (0.025 mol/L):** Opened vial is stable 3 months at 2-8°C.
 4. **Buffer, stock solution (Tris 0.05 mol/L, pH 7.3, 10 mg/L Ciprofloxacin and 1.0% BSA):** Once opened the buffer is stable 3 months at 2-8°C. Prepare a new buffer working solution each day.
 5. **Phospholipid:** Opened vial is stable for 3 months at 2-8°C. Shake gently before use.

Reagents and materials required but not provided
 1. Deionized water, filtered through 0.22 μm or NCCLS type II water.¹¹
 2. Acetic acid 20% or citric acid 2%.
 3. Control Plasma Abnormal and Normal calibrated against an International Standard for Factor VIII
 4. Calibration plasma calibrated against an International Standard
 5. Photometer, 405 nm (and 490 nm for microplate procedure)
 6. Heat incubator 37°C ± 0.2°C
 7. Semi-micro cuvettes
 8. Centrifuge, 2000xg
 9. Plastic test tubes
 10. Stopwatch
 11. Vortex mixer
 12. Calibrated pipettes

Standard	Plasma	Buffer working-solution	Predilution	Final dilution	Buffer working-solution
%	μL	μL	μL	μL	μL
150	-	undiluted	25	2000	
120	200	50	25	2000	
100	100	50	25	2000	
75	100	100	25	2000	
50	100	200	25	2000	
21	100	600	25	2000	


The assigned percentage values of the standard dilutions are those obtained from a normal plasma containing 1.0 IU factor VIII/mL. In case the factor VIII content of the normal plasma differs from this value, an appropriate correction factor should be used.
Preparation of plasma sample
 Use plastic test tubes.
 Test plasma or concentrate 25 μL
 Buffer working solution (2-8°C) 3000 μL
 Mix well. Keep at 2-8°C.
 The assay must be performed within 30 minutes after dilution because of the lability of factor VIII.
Assay
 NOTE: The assay should be performed in plastic material.

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Preparation of the Reagents

The reagents are reconstituted according to the specific instrument application. For microplate and test tube techniques:
 1. **S-2765 + I-2581:** Mix 12.0 ml of sterile water with 12.0 mL of sterile water or NCCLS type II water¹¹ to obtain a concentration of 1.8 mmol/l to be reconstituted.
 2. **Factor IXa + Factor X:** Mix 10.0 ml of sterile water with 10.0 mL of sterile water or NCCLS type II water¹¹ to obtain a concentration of 1.8 mmol/l to be reconstituted.
 3. **CaCl₂:** Ready to use.
 4. **Buffer, stock solution:** Dilute 1:10 (1+9) with sterile water with 10.0 mL of sterile water or NCCLS type II water¹¹ to obtain a concentration of 1.8 mmol/l to be reconstituted.
 5. **Phospholipid:** Ready to use.

Reagents and Stability of the Reagents
 The reagents are reconstituted according to the specific instrument application. For microplate and test tube techniques:
 1. **S-2765 + I-2581:** Stability after reconstitution: 3 months at 2-8°C.
 2. **Factor IXa + Factor X:** Stability after reconstitution: 12 hours at 2-8°C. The solution can be stored frozen in aliquots at -20°C (or at lower temperature) for 3 months. Do not refreeze.
 3. **CaCl₂ (0.025 mol/L):** The reagent is stable for 3 months after opening the vial at 2-8°C.
 4. **Buffer, stock solution (Tris 0.05 mol/L, pH 7.3, 10 mg/L Ciprofloxacin and 1.0 % BSA):** After opening the vial, the buffer is stable for 3 months at 2-8°C. Prepare a new buffer working solution each day.
 5. **Phospholipid:** The reagent is stable for 3 months after opening the vial at 2-8°C. Shake gently before use.

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Specimen collection

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate.
 Centrifugation: 2000 x g for 10-20 minutes at 20-25°C. Refer to NCCLS document H21-A4 for further instructions on specimen collection, handling and storage.¹²

Quality Control

Normal and abnormal controls for plasma or concentrates are recommended for reliable quality control.¹³ Assigned values of Controls should be traceable to the International Standard. Periodically within each run a control should be analyzed. The control material should be treated in the same way as a test sample. A range of allowable variation should be established for controls in each laboratory. If a value outside the established control range is obtained, a complete check of calibration, reagents and instrument performance should be made.

Results

Factor VIII results are reported in % activity (100% factor VIII activity is equivalent to 1.0 IU/mL).

Expected values

Range: 49 - 126 % (2 SD, n=121) in a normal healthy population evaluated with Coatest SP Factor VIII (test tube method).

Due to the many variables that may affect results, each laboratory should establish its own normal range, avoiding inadvertent losses of factor VIII activity.

Procedure

NOTE: All conditions included in this package insert refer to the manual method. Detailed settings for the ACL 8000/9000/10000 including instructions for preparation of the reagents are available on request from Chromogenix.

Two ranges of factor VIII are defined (20-150% and 1-20%).
Range 20-150%:
 Prepare a solution of phospholipid-factor IXa and factor X reagent by mixing:
 • 1 volume of phospholipid
 • 5 volumes of factor IXa+factor X reagent
 Keep on ice or at 2-8°C.
 Shake gently just before use.

Calibration
 A standard curve is required for each Coatest SP Factor VIII kit. Normal human plasma, calibrated against an International Standard, is used for preparation of standard dilutions in plastic tubes using pre-cooled buffer working solution according to the table below:

Standard	Plasma	Buffer working-solution	Predilution	Final dilution	Buffer working-solution
%	μL	μL	μL	μL	μL
150	-	undiluted	25	2000	
120	200	50	25	2000	
100	100	50	25	2000	
75	100	100	25	2000	
50	100	200	25	2000	
21	100	600	25	2000	

The assigned percentage values of the standard dilutions are those obtained from a normal plasma containing 1.0 IU factor VIII/mL. In case the factor VIII content of the normal plasma differs from this value, an appropriate correction factor should be used.

Preparation of plasma sample
 Use plastic test tubes.
 Test plasma or concentrate 25 μL
 Buffer working solution (2-8°C) 3000 μL
 Mix well. Keep at 2-8°C.
 The assay must be performed within 30 minutes after dilution because of the lability of factor VIII.
Assay
 NOTE: The assay should be performed in plastic material.

% FVIII	Acid-stopped method	Initial rate method
Phospholipid+ FIXa + FX (2-8°C)	200 μL	200 μL
Test plasma or standard dilution (2-8°C)	100 μL	100 μL
Mix and incubate at 37°C 4-5 min		
CaCl ₂ (37°C)	100 μL	100 μL
Mix and incubate at 37°C exactly 5 min		
S-2765+I-2581 (37°C)	200 μL	200 μL
Mix and incubate at 37°C exactly 5 min		
Acetic acid 20% or citric acid 2% (20-25°C)	100 μL	

Acid-stopped method: Read the absorbance of the sample against a reagent blank (buffer working solution instead of sample) within 4 hours.
 Because of the large dilution of the plasma, no sample blanks have to be included.

Initial rate method: Transfer immediately to a 1 cm semi-micro cuvette (pre-heated to 37°C) and measure the absorbance change at 405 nm.

Range 1-20%:
 Prepare a solution of phospholipid-factor IXa and factor X reagent by mixing:
 • 1 volume of phospholipid
 • 5 volumes of factor IXa+factor X reagent
 Keep on ice or at 2-8°C.
 Shake gently just before use.
 The mixture is stable for 4 hours at 2-8°C or 12 hours on ice.
Calibration
 A standard curve is required for each Coatest SP Factor VIII kit. Normal human plasma, calibrated against the

Standard	Plasma	Puffer-Arbeitslösung	Verdünnung	Endverdünnung	Puffer-Arbeitslösung
%	μL	μL	μL	μL	μL
150	-	undiluted	25	2000	
120	200	50	25	2000	
100	100	50	25	2000	
75	100	100	25	2000	
50	100	200	25	2000	
21	100	600	25	2000	

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